

AIR/WATER SEPARATORS

ALL PLASTIC - AUTO DRAIN



1H20COM



BVD12946



BVD17373

	CAPACITY (GALLONS)	DIMENSIONS (INCHES)	INLET (INCHES)	OUTLET (INCHES)	VOLTAGE	DRAIN TYPE	HIGH LEVEL SWITCH	WEIGHT (LBS)
BVD17373 Order # 3000102	2.5	22.5 x 11.5	0.5	0.5	N/A	Gravity	✓	22
BVD12946 Order # 3000101	12	28.75 x 16.5	1.5	1.5	24V	Gravity	✓	40
1H20COM Order # 3000001	12	33.75 x 16.5	1.5	1.5	110V	Pump	✓	60

Manufactured By:



R.E Morrison Equipment Inc
Mississauga, ON, Canada
1-800-668-8736
www.remequip.com

Base
DENTAL SYSTEMS **VAC**[™]

www.basevacdental.com

Distributed By:



1 800 368 8106
www.sableindustriesinc.com
ISO 13485..2003 Certified
Health Canada Registeres