

DRY SUCTION C-VAC SERIES

TWO OPERATORY COMPACT DENTAL SUCTION



C-VAC 4.10

POWERFUL - COMPACT - QUIET

	MAX # OF (HVE) USERS	SUCTION POWER ("HG)	POWER (HP)	POWER (AMPI)	VOLTAGE	NOISE LEVEL (DBA)	DIMENSIONS	WEIGHT (LBS)
C-VAC 4.10 Order # 2800014	2	25	.44	7.6	110	64	34 x 16.5	84

- NO OIL, NO WATER, CLEAN DRY EXHAUST
- Integrated high efficiency air water separator automatically drains when the vacuum pump is turned off.
 - Inlet, exhaust and drain all positioned on one side for 'wall tight' installation.
 - Easy clean all-plastic tank makes maintenance simple.
 - Designed to work on pipe sizes from 1/2" to 3".

Manufactured By:



R.E Morrison Equipment Inc
Mississauga, ON, Canada
1-800-668-8736
www.remequip.com
ISO 9001:2008 Certified

Base
DENTAL SYSTEMS **VAC**TM

www.basevacdental.com

Distributed By:



1 800 368 8106
www.sableindustriesinc.com
ISO 13485:2003 Certified
Health Canada Registered