

DRY SUCTION HD SERIES

OIL FREE - WATER FREE - 100% DRY



1HD 4.8
Mini Stationary

1HD 4.16, 4.25, 4.40
Standard Stationary

2HD 4.25, 4.40
Duplex Stationary

PORTA 4.8
Portable System

	MAX # OF (HVE) USERS	SUCTION POWER ("HG)	POWER (HP)	POWER (AMP)	VOLTAGE	NOISE LEVEL (DBA)	DIMENSIONS (INCHES)	WEIGHT (LBS)
PORTA 4.8 Order # 2800008	1	25	0.56	4.8	115V, 1 Phase	65	25 x 20 x 12	65
1HD 4.8 Order # 2800007	1	25	0.56	4.8	115V, 1 Phase	69	13 x 12.5 x 14.5	57
1HD 4.16 Order # 2800103	2	25	0.94	5.2	220V, 1 Phase	67	19 x 18 x 26	87
1HD 4.25 Order # 2800101	4	25	1.2	8.8	220V, 1 Phase	65	19 x 18 x 26	94
1HD 4.40 Order # 2800102	6	25	2.0	9.5	220V, 1 Phase	72	19 x 18 x 26	121
2HD 4.25 Order # 2800012	8	25	2.4	2 x 8.8	220V, 1 Phase	74	38 x 18 x 26	2 x 94
2HD 4.40 Order # 2800006	12	25	4.0	2 x 9.5	220V, 1 Phase	74	38 x 18 x 26	2 x 121

Manufactured By:



R.E Morrison Equipment Inc
Mississauga, ON, Canada
1-800-668-8736
www.remequip.com

Base
DENTAL SYSTEMS **VAC**[™]

www.basevacdental.com

Distributed By:



1 800 368 8106
www.sableindustriesinc.com
ISO 13485..2003 Certified
Health Canada Registeres