

# DRY SUCTION HD SERIES

OIL FREE - WATER FREE - 100% DRY



	MAX # OF (HVE) USERS	SUCTION POWER ("HG)	POWER (HP)	POWER (AMP)	VOLTAGE	NOISE LEVEL (DBA)	DIMENSIONS (INCHES)	WEIGHT (LBS)
<b>PORTA 4.8</b> Order # 2800008	1	25	0.56	4.8	115V, 1 Phase	65	25 x 20 x 12	65
<b>1HD 4.8</b> Order # 2800007	1	25	0.56	4.8	115V, 1 Phase	69	13 x 12.5 x 14.5	57
<b>1HD 4.16</b> Order # 2800103	2	25	0.94	5.2	220V, 1 Phase	67	19 x 18 x 26	87
<b>1HD 4.25</b> Order # 2800101	4	25	1.2	8.8	220V, 1 Phase	65	19 x 18 x 26	94
<b>1HD 4.40</b> Order # 2800102	6	25	2.0	9.5	220V, 1 Phase	72	19 x 18 x 26	121
<b>2HD 4.25</b> Order # 2800012	8	25	2.4	2 x 8.8	220V, 1 Phase	74	38 x 18 x 26	2 x 94
<b>2HD 4.40</b> Order # 2800006	12	25	4.0	2 x 9.5	220V, 1 Phase	74	38 x 18 x 26	2 x 121

Manufactured By:



R.E Morrison Equipment Inc  
Mississauga, ON, Canada  
1-800-668-8736  
[www.remequip.com](http://www.remequip.com)

**Base**  
DENTAL SYSTEMS **VAC**<sup>™</sup>

[www.basevacdental.com](http://www.basevacdental.com)

Distributed By:



1 800 368 8106  
[www.sableindustriesinc.com](http://www.sableindustriesinc.com)  
ISO 13485..2003 Certified  
Health Canada Registeres