

# Experts In Dental Suction

BaseVac Vacuum Systems deliver unrivaled suction performance for the most challenging vacuum applications.



## Upgrading from a Wet Vac to a Dry Vac?

### Better Go BaseVac



[www.basevacdental.com](http://www.basevacdental.com)



905.828.6301



[info@basevacdental.com](mailto:info@basevacdental.com)

#### **100% Oil & Water Free**

Save on water costs. No oil changes or oily exhaust residue.

#### **Industry Leading Suction Performance**

Adjustable vacuum up to 25" Hg

#### **Works with 1" Dia or Smaller Pipes**

Replace wet vacuums without replacing in-floor plumbing.

#### **5-Year Warranty**

Unlimited hour hassle-free warranty.

# Dry Vacuum Systems



100% Oil & Water Free | Industry Leading Suction 10”Hg-25”Hg | Rotary Vane Technology



Model #	2M4.16	2M4.25	1HD4.25	1HD4.40	2HD4.25	2HD4.40	3HD4.40	4HD4.40
Max Users	2	2.5	2.5	4	5	8	12	16
Suction Power	10-25 "Hg	10-25 "Hg	10-25 "Hg	10-25 "Hg	10-25 "Hg	10-25 "Hg	10-25 "Hg	10-25 "Hg
Power (Hp)	0.94 Hp	1.2 Hp	1.2 Hp	2 Hp	3.4 Hp	4 Hp	6 Hp	8 Hp
FLA (Amps)	5.2 A	8.8 A	6.9A	10.5A	6.9A	10.5A	10.5A	10.5A
Voltage Rating (Min - Max)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)	230V, 1 Phase (207-253V)
Noise Level	67 dBA	67 dBA	67 dBA	69 dBA	70 dBA	72 dBA	75 dBA	78 dBA
Dimensions	L 17" x W 17" x H 30"	L 17" x W 17" x H 39"	L 28" x W 18" x H 19"	L 28" x W 18" x H 19"	L 28" x W 18" x H 38"	L 28" x W 18" x H 38"	L 28" x W 18" x H 57"	L 28" x W 18" x H 76"
Weight (lbs)	86 lbs	98 lbs	94 lbs	121 lbs	188 lbs	242 lbs	363 lbs	484 lbs
Rec. Vacuum Lines	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"
Min. Vacuum Lines	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"
# of Motors	1	1	1	1	2	2	3	4
Upgradeable	No	No	Yes	Yes	Yes	Yes	Yes	Yes
Buck-Boosts (Optional)	1	1	1	1	2	2	3	4

# Oil-Free Compressors

Operating PSI: 70-100 | Oil-less Air Compressors | PERMA-DRY Air Purification | Quiet Operation



Model #	S1.5	S2.0	S1.5E	S2.0E	D1.5	D2.0
Users	3	4-5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	6-7	8-10
Operatories	3-4	3-5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	5-7	8-10
CFM @ 80PSI	5	6.5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	10	14
Power (hp)	1.5 hp	2 hp	Same as Non-Enclosure Model	Same as Non-Enclosure Model	2 x 1.5 hp	2 x 2 hp
Power (Amps)	8.5 A	9 A	Same as Non-Enclosure Model	Same as Non-Enclosure Model	17 A	18 A
Electrical	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz
Dimensions	L 26" x W 20" x H 25"	L 26" x W 25" x H 30"	L 29" x W 25" x H 32"	L 29" x W 25" x H 32"	L 35" x W 24" x H 36"	L 35" x W 24" x H 36"
Noise Level	72 dBA	74 dBA	54 dBA	54 dBA	74 dBA	74 dBA
Weight (lbs)	132 lbs	154 lbs	172 lbs	172 lbs	231 lbs	242 lbs
Tank	13 Gal.	13 Gal.	13 Gal.	13 Gal.	20 Gal.	20 Gal.
Silent Option (-20 dBA)	Yes	Yes	N/A	N/A	No	No

# Equipment Room Packages

Easy Installation | Industry Leading Suction | Oil and Water Free



Dry Vacuum System



Oil-Free Air Compressor



Buck-Boost Transformer



Remote Switch

## 2+ User Package

1HD4.25 Vacuum | 2.5 User | #2800004  
S1.5 Compressor | 3 User | #2880106  
Dual Remote Switch | #2900103  
0.25 KVA Buck Boost (x2) | #2900601

## 4 User Package

1HD4.40 Vacuum | 4 User | #2800005  
S2.0 Compressor | 4-5 User | #2880109  
Dual Remote Switch | #2900103  
0.25 KVA Buck Boost | #2900601  
0.5 KVA Buck Boost | #2900602

## 6+ User Package

2HD4.40 Vacuum | 8 User | #2800006  
D1.5 Compressor | 6-7 User | #2880116  
Triple Remote Switch | #2900106  
0.5 KVA Buck Boost (x3) | #2900602

## 8 User Package

2HD4.40 Vacuum | 8 User | #2800006  
D2.0 Compressor | 8-10 User | #2880111  
Triple Remote Switch | #2900106  
0.5 KVA Buck Boost (x3) | #2900602

## 10+ User Package

3HD4.40 Vacuum | 12 User | #2800030  
D2.0 Compressor | 8-10 User | #2880111  
Quad Remote Switch | #2880199  
0.5 KVA Buck Boost (x4) | #2900602



# BaseVac Testimonials

BaseVac is trusted by top professionals across North America



I can highly recommend this Canadian made system to all of our colleagues, as the BaseVac system is of high quality and will exceed the needs of any dental clinic.

## Dr. Kate Bazydlo

Winston Churchill Dental  
Mississauga, ON



## Dr. Bob Lail

Dr. Lail Associates  
Woodstock, ON

There is no question; the quality, service, and support from a BaseVac is why I would recommend their systems to any Dental Professional looking for mechanical room equipment.



## Dr. D. Keleman

King East Dental  
Stoney Creek, ON

Not only is their product offering superior to other vacuums on the market but, their sales staff is extremely helpful in ensuring that I had the system that met my unique application.

**BaseVac Dental products are sold exclusively by  
Professional Sales Associates (PSA)**



Professional Sales Associates Inc. (PSA)

5045 Park Ave West, Ste. 1B

Seville, OH | 44273 | U.S.A

330-299-7343 | [info@profsales.com](mailto:info@profsales.com) | [www.profsales.com](http://www.profsales.com)