

Experts In Dental Suction

BaseVac Vacuum Systems deliver unrivaled suction performance for the most challenging vacuum applications.



Upgrading from a Wet Vac to a Dry Vac?

Better Go BaseVac



www.basevacdental.com



905.828.6301



info@basevacdental.com

100% Oil & Water Free

Save on water costs. No oil changes or oily exhaust residue.

Industry Leading Suction Performance

Adjustable vacuum up to 25" Hg

Works with 1" Dia or Smaller Pipes

Replace wet vacuums without replacing in-floor plumbing.

5-Year Warranty

Unlimited hour hassle-free warranty.

Dry Vacuum Systems



100% Oil & Water Free | Industry Leading Suction 10"Hg-25"Hg | Rotary Vane Technology









Model #	2M4.16	2M4.25	1HD4.25	1HD4.40	2HD4.25	2HD4.40	3HD4.40	4HD4.40
Max Users	2	2.5	2.5	4	5	8	12	16
Suction Power	10-25 "Hg							
Power (Hp)	0.94 Hp	1.2 Hp	1.2 Hp	2 Hp	3.4 Hp	4 Hp	6 Hp	8 Hp
FLA (Amps)	5.2 A	8.8 A	6.9A	10.5A	6.9A	10.5A	10.5A	10.5A
Voltage Rating (Min - Max)	230V, 1 Phase (207-253V)							
Noise Level	67 dBA	67 dBA	67 dBA	69 dBA	70 dBA	72 dBA	75 dBA	78 dBA
Dimensions	L 17" x W 17" x H 30"	L 17" x W 17" x H 39"	L 28" x W 18" x H 19"	L 28" x W 18" x H 19"	L 28" x W 18" x H 38"	L 28" x W 18" x H 38"	L 28" x W 18" x H 57"	L 28" x W 18" x H 76"
Weight (lbs)	86 lbs	98 lbs	94 lbs	121 lbs	188 lbs	242 lbs	363 lbs	484 lbs
Rec. Vacuum Lines	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"	1 1/2"
Min. Vacuum Lines	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"
# of Motors	1	1	1	1	2	2	3	4
Upgradeable	No	No	Yes	Yes	Yes	Yes	Yes	Yes
Buck-Boosts (Optional)	1	1	1	1	2	2	3	4

Oil-Free Compressors

Operating PSI: 70-100 | Oil-less Air Compressors | PERMA-DRY Air Purification | Quiet Operation







Model #	S1.5	\$2.0	S1.5E	\$2.0E	D1.5	D2.0
Users	3	4-5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	6-7	8-10
Operatories	3-4	3-5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	5-7	8-10
CFM @ 80PSI	5	6.5	Same as Non-Enclosure Model	Same as Non-Enclosure Model	10	14
Power (hp)	1.5 hp	2 hp	Same as Non-Enclosure Model	Same as Non-Enclosure Model	2 x 1.5 hp	2 x 2 hp
Power (Amps)	8.5 A	9 A	Same as Non-Enclosure Model	Same as Non-Enclosure Model	17 A	18 A
Electrical	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz	220V / 1 / 60Hz
Dimensions	L 26" x W 20" x H 25"	L 26" x W 25" x H 30"	L 29" x W 25" x H 32"	L 29" x W 25" x H 32"	L 35" x W 24" x H 36"	L 35" x W 24" x H 36"
Noise Level	72 dBA	74 dBA	54 dBA	54 dBA	74 dBA	74 dBA
Weight (lbs)	132 lbs	154 lbs	172 lbs	172 lbs	231 lbs	242 lbs
Tank	13 Gal.	13 Gal.	13 Gal.	13 Gal.	20 Gal.	20 Gal.
Silent Option (-20 dBA)	Yes	Yes	N/A	N/A	No	No

Equipment Room Packages



Easy Installation | Industry Leading Suction | Oil and Water Free







Oil-Free Air Compressor



Remote Switch

2+ User Package

1HD4.25 Vacuum | 2.5 User | #2800004 S1.5 Compressor | 3 User | #2880106 Dual Remote Switch | #2900103 0.25 KVA Buck Boost (x2) | #2900601

6+ User Package

2HD4.40 Vacuum | 8 User | #2800006 D1.5 Compressor | 6-7 User | #2880116 Triple Remote Switch | #2900106 0.5 KVA Buck Boost (x3) | #2900602

4 User Package

1HD4.40 Vacuum | 4 User | #2800005 S2.0 Compressor | 4-5 User | #2880109 Dual Remote Switch | #2900103 0.25 KVA Buck Boost | #2900601 0.5 KVA Buck Boost | #2900602

8 User Package

2HD4.40 Vacuum | 8 User | #2800006 D2.0 Compressor | 8-10 User | #2880111 Triple Remote Switch | #2900106 0.5 KVA Buck Boost (x3) | #2900602

10+ User Package

3HD4.40 Vacuum | 12 User | #2800030 D2.0 Compressor | 8-10 User | #2880111 Quad Remote Switch | #2880199 0.5 KVA Buck Boost (x4) | #2900602

Base VAC

BaseVac Testimonials

BaseVac is trusted by top professionals across North America



BaseVac offered a unique design for my specific application. Now thanks to their team, our suction is stronger, more consistent, and more efficient than it has ever been before."

Dr. Anna FromzelStamford Family Dental Stamford, CT





mechanical room."

Dr. Ali MiameeReston Family Dental
Reston, VA

The four years of powerful vacuum and air we received is why I only trust BeaseVac Dental systems in my



Honestly they are great. Up to this point, I have had zero issues. They work. They do the job. We've had alternatives but so far, BaseVac has outperformed all of those."

BaseVac Dental products are sold exclusively by Professional Sales Associates (PSA)



Professional Sales Associates Inc. (PSA)
5045 Park Ave West, Ste. 1B
Seville, OH | 44273 | U.S.A
330-299-7343 | info@profsales.com | www.profsales.com